



## Simplify Patient Positioning

Positioning patients for radiography examinations is convenient, efficient and precise with the mobile CARESTREAM SPT1 Radiographic Table. Designed for maximum maneuverability, this versatile imaging table is suitable for DR, CR and film-based procedures.

## ADVANTAGES AT A GLANCE

- Multi-functional table with 4-way floating top
- Lightweight and easy to move
- Pedestal design allows the table to be freely positioned at any angle
- Easy access float and wheel brakes



# CARESTREAM SPT1 Radiographic Table



## Easy to Handle

- Directional wheel castors and a large handgrip allow for easy movement in any direction.
- Hydraulic control elevates the table with no external power source or cable.

## Easy to Clean

- Streamlined table design simplifies cleaning and disinfecting.

## Maintenance-Free

- This rugged, self-contained table requires no maintenance.

## Floating Tabletop

- A four-way floating tabletop moves easily for precise patient positioning.
- Foot-pedal control locks the tabletop in position.
- Tabletop width enhances positioning ease.

## Four-Wheel Brakes

- A brake on every wheel helps to ensure stability and patient safety.
- All brakes are set and released simultaneously with one press of a foot pedal. Three pairs of pedals are conveniently located.
- The table can be anchored to the floor for increased collision avoidance.

## Product Specifications

### • Table Dimensions

Length (including handle): 85.8 in. (218 cm)  
Width: 28.7 in. (73 cm)  
Maximum height: 33.9 in. (86 cm)

### • Tabletop Dimensions

Length: 83.5 in. (212 cm)  
Width: 25.6 in. (65 cm)  
Thickness: 2.1 in. (5.4 cm)

### • Threshold Height

Maximum Height: 0.6 in. (15 mm)

### • Table Weight

341.7 lbs (155 kg)

### • Allowable Patient Weight

Maximum patient weight: 500 lbs (225 kg)

### • X-ray Transparent Surface

Length: 83.5 in. (212 cm)  
Width: 25.6 in. (65 cm)

### • X-ray Absorption Factor

Factor:  $\leq 0.8$  mm aluminum equivalence

### • Travel Position of the Tabletop

Elevating range: 9.8 in. (25 cm)  
Longitudinal float: 22.5 in. (56 cm)  
Transversal float:  $\pm 4.8$  in. ( $\pm 12$  cm)

## Options

- Pivot Kit provides the components needed to allow you to rotate the table around a defined point in the floor.
- Table Pad is available to reduce patient discomfort.
- Table Compression Band enables quality imaging by reducing thickness and/or restricting movement of patients.

## Environmental Requirements

### • Transit and Storage

The receiving and storage area(s) must be dry and able to provide the proper humidity and temperature control required for the equipment.

### • Storage Temperature

-9 to 150°F (-23 to 66°C)

### • Relative Humidity

35 to 75% (non-condensing)

### • Atmospheric Pressure

70 to 106 kPa

### • Operating Temperature

50 to 104°F (10 to 40°C)

